

IHE Integration Statement			
Vendor	Product Name	Version	Date
InterOperability Bidco, Inc.	Corepoint Integration Engine	7.3	1/26/18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Node	None
Audit Trail and Node Authentication	Secure Node	FQDN Validation of Server Certificate
Audit Trail and Node Authentication	Secure Node	ATX: FHIR FEED
Audit Trail and Node Authentication	Secure Node	STX:TLS 1.2 Floor using BCP195
Consistent Time	Time Client	None
Mobile Care Services Discovery	Care Services Selective Consumer	None
Mobile access to Health Documents	Document Consumer	None
Mobile access to Health Documents	Document Recipient	None
Mobile access to Health Documents	Document Source	None
Patient Demographics Query for Mobile	Patient Demographics Consumer	None
Patient Demographics Query for Mobile	Patient Demographics Supplier	None
Patient Location Tracking	Patient Location Tracking Consumer	None
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Simple Observations
Query for Existing Data for Mobile	Clinical Data Consumer	None
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Allergies and Intolerances

Internet address for vendor's IHE information :
<https://corepointhealth.com/company/corepoint-health-ihe-integration-statement/>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
InterOperability Bidco, Inc.	Corepoint Integration Engine	7.3	1/26/18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Conditions
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Diagnostic Results
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Medications
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Immunizations
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Procedures
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Encounters
Query for Existing Data for Mobile	Clinical Data Consumer	QEDm Document Provenance
Retrospective Data Query	Retrospective Data Consumer	None
Remove Metadata and Documents	Document Administrator	Remote Registry
Remove Metadata and Documents	Document Administrator	None
Remove Metadata and Documents	Document Administrator	Remote Repository
Restricted Metadata Update	Update Initiator	None
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Source	None

Internet address for vendor's IHE information :
<https://corepointhealth.com/company/corepoint-health-ihe-integration-statement/>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net