

IHE Integration Statement			
Vendor	Product Name	Version	Date
InterOperability Bidco, Inc.	Rhapsody Integration Engine	6.5	10/30/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Node	STX:TLS 1.2 Floor using BCP195
Audit Trail and Node Authentication	Secure Node	None
Audit Trail and Node Authentication	Secure Node	ATX: FHIR FEED
Consistent Time	Time Client	None
Device Enterprise Communication	Device Observation Consumer	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query	Patient Demographics Consumer	Patient Demographics and Visit Query
Patient Demographics Query for Mobile	Patient Demographics Consumer	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	PIX Update Notification
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	PIXV3 Update Notification Transaction
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None

Internet address for vendor's IHE information :
<https://corepointhealth.com/company/corepoint-health-ihe-integration-statement/>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
InterOperability Bidco, Inc.	Rhapsody Integration Engine	6.5	10/30/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Community Access	Document Consumer	Asynchronous Web Services Exchange (WS-Addressing based)
Cross-Community Access	Document Consumer	None
Cross-Community Access	Initiating Gateway	AS4 Asynchronous Web Services Exchange
Cross-Community Access	Initiating Gateway	None
Cross-Enterprise Document Reliable Interchange	Document Recipient	None
Cross-Enterprise Document Reliable Interchange	Document Recipient	Accepts Limited Metadata
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Consumer	On-Demand Documents
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Document Source	Folder Management
Cross-Enterprise Document Sharing	Document Source	Document Addendum
Cross-Enterprise Document Sharing	Document Source	Document Transformation

Internet address for vendor's IHE information :
<https://corepointhealth.com/company/corepoint-health-ihe-integration-statement/>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
InterOperability Bidco, Inc.	Rhapsody Integration Engine	6.5	10/30/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None

Internet address for vendor's IHE information :
<https://corepointhealth.com/company/corepoint-health-ihe-integration-statement/>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net